
## **UAB "LT Biotech"**



Į.K. 302303586, PVM kodas LT 100004741118, Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į. Registrų centras Vilniaus filialas Tel/fax +370 5 216 02 27

# Cyropreservation Medium with Fetal Bovine Serum Cat. No.: MFFBS50ML (50 ml)

#### **General Information**

Cryoprotective medium that contains 10 % DMSO and 20 % FBS. The serum used is extensively tested to protect cells during cell preservation.

## **Applications**

- Cryopreservation of a wide range of cell types with high viability
- Ready-to-use solution

Appearance	Clear red liquid
Storage and shelf life	Store at ≤-15°C. is a light sensitive solution. It should be protected from light during storage.
Shipping conditions	Frozen (Dry Ice)
Thawing	+37°C water bath or overnight at +2°C to +8°C. Swirl gently to homogenize.

#### Instructions for Use

#### Freezing Protocols

Before cryopreservation cells should be checked for contamination. Can be used with any standard freezing protocol.


### Cryopreservation of Suspension Cultures

- 1. Count the number of viable cells to be cryopreserved. Cells should be in mid-log phase of growth. Centrifuge the cells for 5 min to pellet cells (200 to 400 g). Remove the supernatant down to the smallest volume without disturbing the cells.
- 2. Resuspend cells in pre-cooled ( $+4^{\circ}$ C to  $+8^{\circ}$ C) to a concentration of  $5x10^{6}$  to  $1x10^{7}$  cells/ml.
- 3. Aliquot into cryogenic storage vials. Place vials at +4°C and start the freezing procedure within 5 min. Cells are frozen slowly at +1°C/min (by programmable coolers or by placing vials in an insulated box in a -70°C to -90°C freezer).
- 4. Then transfer storage vials to liquid nitrogen storage.

#### Cryopreservation of Adherent Cultures

- 1. Detach cells from the substrate with a gentle dissociating agent. Especially with sensitive cells use Accutase\* (Cat. No. ACC-1B) to avoid cell damage. Inactivate dissociating agent if necessary.
- 2. Resuspend the detached cells in complete growth medium and establish the viable cell count.
- 3. Centrifuge for 5 min to pellet cells (200 to 400 g). Remove the supernatant down to the smallest volume without disturbing the cells.
- 4. Resuspend cells in pre-cooled (+4°C to +8°C) to a concentration of 5x106 to 107 cells/ml.

## **UAB "LT Biotech"**



Į.K. 302303586, PVM kodas LT 100004741118, Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į. Registrų centras Vilniaus filialas Tel/fax +370 5 216 02 27

- 5. Aliquot into cryogenic storage vials. Place vials at +4°C and start the freezing procedure within 5 min. Cells are frozen slowly at +1°C/min (by programmable coolers or by placing vials in an insulated box in a -70°C to -90°C freezer).
- 6. Then transfer storage vials to liquid nitrogen storage.

## Thawing of Cryopreserved Cells

Cryopreserved cells can be thawed by the following procedures:

#### Centrifugation

- 1. Remove cells from storage and thaw quickly in a +37°C water bath. Capricorn Scientific recommends eye protection by using approved safety goggles. We also suggest the use of safety gloves to protect uncovered skin.
- 2. Place 1 to 2 ml of thawed cells in ~25 ml of complete growth medium. Mix cell suspension gently.
- 3. Centrifuge the cells at ~80 g for 2 to 3 min.
- 4. Check clarity of the supernatant and visibility of a consolidated cell pellet. Discard supernatant without disturbing the cells.
- 5. Gently resuspend the cells in complete growth medium and perform a viable cell count.
- 6. Plate the cells. Cell inoculum should be at least 3x10<sup>5</sup> viable cells/ml.

#### Direct plating

- 1. Remove cells from storage and thaw quickly in a +37°C water bath. Capricorn Scientific recommends eye protection by using approved safety goggles. We also suggest the use of safety gloves to protect uncovered skin
- 2. Plate cells directly with complete growth medium. Use 10 to 20 ml of complete medium per 1 ml of frozen cells. Cell inoculum should be at least 3x10<sup>5</sup> cells/ml.
- Culture cells for 12 to 24 h. Replace medium with fresh complete growth medium to remove cryopreservative.

We recommend thawing procedure 1, especially when handling sensitive cells.

#### **Precautions and Disclaimer**

This product is for research use only.

**Caution:** contains Dimethyl Sulfoxide (DMSO). Do not breathe gas/fumes/vapour/spray. Avoid contact with eyes and skin. Irritant to eyes, respiratory system and skin. S23 S24/25.