



LT BIOTECH LTD,
Mokslininkų g. 6A,
LT-08412 Vilnius,
Lithuania
www.ltbiotech.lt

LTAf06 Defoamer Specification

LTAf06 silicone defoamer is mainly recommended for food and drug fermentation and processing.

FRATURES

The substance components of the antifoam do not contain the substances listed in the SVHC list of the REACH Regulation. Antifoaming agent component does not contain NP, APEO.

COMPOSITION

Composed of modified polysiloxane, synergist and silica particles.

| No | Items | Specification |
|----|------------------------|--------------------------------------------------------------------------------------------|
| 1 | Appearance | White or yellowish uniform liquid emulsion, without any sediment and mechanical impurities |
| 2 | PH value | 6.5~8.5 |
| 3 | Nonvolatile content, % | 10.0±1.0 |

HOW TO USE

LTAf06 Antifoam can be diluted with thickening water. In the process of dilution, the speed of stirring is between 30 to 60 rpm and the time of stirring is recommended for no more than 10 minutes. It is suggested that the product is added in a place where the antifoam can disperse easily.

PACKAGE AND STORAGE

The product is available in 25kg or 200kg plastic buckets; sizes can be customized according to special requirements. This product can be stored as a normal chemical, away from heat sources and direct sunlight. No exposure to rain. Keep away from flame. Don't add acid, alkali, salt and other corrosive substances to this product. Ensure the container is sealed after using. If this product begins to stratify after long storage, stir it evenly at 30 to 60 rpm. This product may freeze under 0°C; if this happens, stir well after defrosting the antifoam at a rate of 30 to 60 rpm. In these cases, stirring will not affect the performance of the product if carried out at the recommended speed. When stored between 5 and 40°C in the original unopened containers, the product is usable 9 months from the date of manufacture.

LIMITIED WARRANTY INFORMATION

Please read carefully: The information contained herein is offered from experiences and knowledge in our company. The proposed application method of this product is for reference only. It is recommended that the customer conduct preliminary tests, ensuring product safety and effectiveness. Each customer has different system conditions.