

UAB "LT Biotech"

Į.K. 302303586, PVM kodas LT 100004741118,
Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į.
Registru centras Vilniaus filialas
Tel/fax +370 5 216 02 27


DMEM/Ham'sF-12,w/oL-Glutamine,with15mMHEPESLTM-B02

General Information

Appearance	Clearredorangesolution
CO ₂ concentration, optimum	5.5%
Storage and shelf life	Storeat+2°C–+8°Cprotectedfromlight. Onceopened,storeat4°Candusewithin6-8weeks.
Shipping conditions	Ambient

Formulation

Components	Concentration mg/l		
Amino Acids		Vitamin B ₁₂	0.68
Glycine	18.75	Inorganic Salts	
L-Alanine	4.45	CaCl ₂ 2 H ₂ O	154.50
L-Arginine HCl	147.50	CuSO ₄ 5 H ₂ O	0.0013
L-Asparagine H ₂ O	7.50	Fe (NO ₃) ₃ 9 H ₂ O	0.05
L-Aspartic Acid	6.65	FeSO ₄ 7 H ₂ O	0.417
L-Cysteine HCl H ₂ O	17.56	KCl	311.80
L-Cystine 2 HCl	31.29	MgCl ₂ 6 H ₂ O	61.20
L-Glutamic Acid	7.35	MgSO ₄	48.84
L-Histidine HCl H ₂ O	31.48	NaCl	6996.00
L-Isoleucine	54.47	NaHCO ₃	2438.00
L-Leucine	59.05	Na ₂ HPO ₄	71.02
L-Lysine HCl	91.25	NaH ₂ PO ₄	54.30
L-Methionine	17.24	ZnSO ₄ 7 H ₂ O	0.43
L-Phenylalanine	35.48	Other Components	
L-Proline	17.25	D-Glucose	3151.00
L-Serine	26.25	HEPES	3574.50
L-Threonine	53.45	Hypoxanthine	2.10
L-Tryptophan	9.02	Linoleic Acid	0.04
L-Tyrosine 2 Na 2 H ₂ O	55.79	Lipoic Acid	0.105
L-Valine	52.85	Phenol Red Sodium Salt	8.63
Vitamins		Putrescine 2 HCl	0.081
D-Biotin	0.0035	Sodium Pyruvate	55.00
Choline Chloride	8.98	Tymidine	0.365
D-Calcium Pantothenate	2.24		
Folic Acid	2.66		
myo-Inositol	12.60		
Nicotinamide	2.02		
Pyridoxine HCl	0.31		
Riboflavin	0.22		
Thiamine HCl	2.17		



UAB "LT Biotech"

J.K. 302303586, PVM kodas LT 100004741118,
Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į.
Registru centras Vilniaus filialas
Tel/fax +370 5 216 02 27

Precautions and Disclaimer

This product is for research use only.