

UAB "LT Biotech"

Į.K. 302303586, PVM kodas LT 100004741118,
Rugiu 21-24, LT-08419, Vilnius Reg. Nr. 127918,
V.Į. Registrų centras Vilniaus filialas
Tel/fax +370 5 216 02 27


RPMI 1640, w/o L-Glutamine, with 25 mM HEPES LTLM-B25

General Information

Appearance	Clear red orange solution
CO ₂ concentration, optimum	4.5 %
Storage and shelf life	Store at +2°C – +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.
Shipping conditions	Ambient

Formulation

Components	Concentrationmg/l		
		D-Biotin	0.20
Amino Acids		CholineChloride	3.00
Glycine	10.00	D-Calcium Pantothenate	0.25
L-Arginine	200.00	FolicAcid	1.00
L-Asparagine	50.00	myo-Inositol	35.00
L-Aspartic Acid	20.00	Nicotinamide	1.00
L-Cystine 2HCl	65.20	Pyridoxine HCl	1.00
L-Glutamic Acid	20.00	Riboflavin	0.20
L-Histidine	15.00	Thiamine HCl	1.00
L-Hydroxy-L-Proline	20.00	Vitamin B ₁₂	0.005
L-Isoleucine	50.00	Inorganic Salts	
L-Leucine	50.00	Ca(NO ₃) ₂ ·4H ₂ O	100.00
L-Lysine HCl	40.00	KCl	400.00
L-Methionine	15.00	MgSO ₄	48.84
L-Phenylalanine	15.00	NaCl	6000.00
L-Proline	20.00	NaHCO ₃	2000.00
L-Serine	30.00	Na ₂ HPO ₄	800.00
L-Threonine	20.00	Other Components	
L-Tryptophan	5.00	D-Glucose	2000.00
L-Tyrosine 2Na 2H ₂ O	28.83	HEPES	5958.00
L-Valine	20.00	L-Glutathione Reduced	1.00
Vitamins		Phenol Red Sodium Salt	5.30
p-Amino Benzoic Acid	1.00		



UAB "LT Biotech"

Į.K. 302303586, PVM kodas LT 100004741118,
Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918,
V.Į. Registrų centras Vilniaus filialas
Tel/fax +370 5 216 02 27

Precautions and Disclaimer

This product is for research use only.