

UAB "LT Biotech"

Į.K. 302303586, PVM kodas LT 100004741118,
Rugiu 21-24, LT-08419, Vilnius Reg. Nr. 127918,
V.Į. Registrų centras Vilniaus filialas
Tel/fax +370 5 216 02 27


RPMI 1640, w/o L-Glutamine, w/o Phenol Red LTLM-B28

General Information

Appearance	Clear red orange solution
CO ₂ concentration, optimum	4.5 %
Storage and shelf life	Store at +2°C – +8°C protected from light. Once opened, store at 4° C and use within 6-8 weeks.
Shipping conditions	Ambient

Formulation

Components	Concentration mg/l	Vitamins	
Amino Acids		p-Amino Benzoic Acid	1.00
Glycine	10.00	D-Biotin	0.20
L-Arginine	200.00	Choline Chloride	3.00
L-Asparagine H ₂ O	56.81	D-Calcium Pantothenate	0.25
L-Aspartic Acid	20.00	Folic Acid	1.00
L-Cystine 2HCl	65.20	myo-Inositol	35.00
L-Glutamic Acid	20.00	Nicotinamide	1.00
L-Histidine	15.00	Pyridoxine HCl	1.00
L-Hydroxy-L-Proline	20.00	Riboflavin	0.20
L-Isoleucine	50.00	Thiamine HCl	1.00
L-Leucine	50.00	Vitamin B ₁₂	0.005
L-Lysine HCl	40.00	Inorganic Salts	
L-Methionine	15.00	Ca(NO ₃) ₂ ·4 H ₂ O	100.00
L-Phenylalanine	15.00	KCl	400.00
L-Proline	20.00	MgSO ₄	48.84
L-Serine	30.00	NaCl	6000.00
L-Threonine	20.00	NaHCO ₃	2000.00
L-Tryptophan	5.00	Na ₂ HPO ₄	800.00
L-Tyrosine 2Na 2H ₂ O	28.83	Other Components	
L-Valine	20.00	D-Glucose	2000.00
		L-Glutathione Reduced	1.00



UAB “LT Biotech”

Į.K. 302303586, PVM kodas LT 100004741118,
Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918,
V.Į. Registrų centras Vilniaus filialas
Tel/fax +370 5 216 02 27

Precautions and Disclaimer

This product is for research use only.