



UAB "LT Biotech"

J.K. 302303586, PVM kodas LT 100004741118,
Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918,
V.Į. Registrų centras Vilniaus filialas
Tel/fax +370 5 216 02 27

Adult Bovine Serum #LTABS100 #LTABS500#TABS1000

General Information

Adult Bovine Serum is obtained from animals which are older than 12 months and contains more immunoglobulins than Fetal Bovine Serum and has increased protein content. Adult Bovine Serum is a cost effective alternative to Fetal Bovine Serum.

Upon request this serum can be offered with treatments such as heat inactivation or gamma irradiation.

| | |
|------------------------|---|
| Appearance | Clear amber-brownish liquid |
| Storage and shelf life | Store at $\leq -15^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened, store at 4°C and use within 4-6 weeks. |
| Shipping conditions | Frozen (Dry ice) |
| Thawing | $+37^{\circ}\text{C}$ water bath or overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$. Swirl gently to homogenize. |

Quality Control

Only serum batches which pass our strict quality control are released for sale. Standard parameters which are determined include pH, osmolality, content of protein, albumin, hemoglobin, endotoxin level, sterility, mycoplasma detection and virus testing.

Precautions and Disclaimer

This product is for research use only.