


## Datasheet

<b>Catalogue #</b>	LT-AG-SPA-01
<b>Product name</b>	Staphylococcal protein A (SPA)
<b>Source</b>	from recombinant <i>E. coli</i> strain
<b>Specificity</b>	React specifically with Fc region of IgG some mammals species
<b>Molecular weight</b>	33,0 kD
<b>Theoretical pI</b>	6,06
<b>Purification tag</b>	6xHis
<b>Unit</b>	1 mg, 5 mg, 10 mg
<b>Concentration</b>	Lot dependent
<b>Purification</b>	chromatography
<b>Buffer</b>	PBS
<b>Stabilizer</b>	none, may be added on request
<b>Preservative</b>	None
<b>Application</b>	ELISA, Lateral Flow
<b>Sterility</b>	Filter sterilized
<b>Long Term Storage</b>	Lyophilized
<b>Material Safety note</b>	This product is sold as an antigen preparation for research use only. Standard Laboratory Practices should be followed when handling this material.
<b>Validation data:</b>	

SDS-PAGE



Coomassie stained SPA preparation after resolution in 12% PAGE at 35 mA 1h