

Į.K. 302303586, PVM kodas LT 100004741118, Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į. Registrų centras Vilniaus filialas Tel/fax +370 5 216 02 27


CHO CHO FIRST Media System[™] Cat. No.: Starter Kit (LTCHO start); Expansion Kit (LTCHO start 2)

The CHO First [™] consists of two media and two feed supplements that are built out of the same basis of raw mate-rials. All parts of the system are animal component free and do not contain any hydrolysates. Thus, the CHO First[™] is fully chemically defined.


- CHO First A for selective adaption not rich in nutrients but balanced to adapt efficiently
- CHO First E for production non selective and rich in nutrients for an efficient production
- CHO First Feed 1 & CHO First Feed 2 supplements rich in nutrients for an balanced predetermined feeding strategy

This innovative combination allows an easy installation, easy feeding & easy scale-up process for achieving the best results.

Best viability with CHO First ™



The media and supplements of the CHO First [™] contain nutrients that recombinant CHO cells need for good growth and a higher yield of recombinant proteins with high potency. Fed-batch process was stopped when viability dropped below 80%.



Highest production yield with CHO First ™



Į.K. 302303586, PVM kodas LT 100004741118, Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į. Registrų centras Vilniaus filialas Tel/fax +370 5 216 02 27

The excellent performance is based on the unique Media System with balanced media and supplements. Thus, it is possible to scale-up the production process without any adaptation. Due to the extreme healthy conditions the highest production yield can be reached.

CHO First™ Components

Product	Volume	Cat. No.
CHO First ™, Starter Kit		LTCHO start
CHO First ™, Expansion Kit		LTCHO start 2
CHO First A	500 ml	LTCHOF/A500
CHO First E	500 ml / 1000 ml	LTCHOF/E500 / CLTCHOF/E1000
CHO First Feed 1	100 ml / 500 ml	LTCHOF/F1.100/ LTCHOF/F1.500
CHO First Feed 2	10 ml / 50 ml	CLTCHOF/F2.10 / CLTCHOF/F2.5

Test Kit Composition

Liquid Kit Component	Starter Kit	Expansion Kit
CHO First A	1 x 500 ml	-
CHO First E	1 x 500 ml	5 x 1000 ml
CHO First Feed 1	2 x 100 ml	5 x 500 ml
CHO First Feed 2	2 x 10 ml	5 x 50 ml

Related Products

Product	Volume	Cat. No.
L-Glutamine (200 mM)	100 ml	LTGA02
Recombinant Insulin (5 mg/ml)	5 ml	HRINS5ML
Glucose (250 g/l)	50 ml	GLUS50ML