
UAB "LT Biotech"



Į.K. 302303586, PVM kodas LT 100004741118, Rugių 21-24, LT-08419, Vilnius Reg. Nr. 127918, V.Į. Registrų centras Vilniaus filialas Tel/fax +370 5 216 02 27

Hanks'BalancedSalts,w/oCa&Mg,withPhenolRed,sterile-filtered LTBSS03

General Information

Hank'sBalancedSaltSolution(HBSS)isusedforavarietyofcellcultureapplications, such as rinsingcells, transportingcells ortissues amples, dilutingcells for counting, and preparing reagents. The essential function of abalanced saltsolution is to maintain pH and os motic balance as well as providing cells with water and essential inorganic ions.

Product Specifications

Appearance	Clear red liquid	
Storage and shelf life	Storeat+2to+8°Cuponarrival.	
	Onceopenedstoreat+2to+8°Candusewithin24month.	
Shipping conditions	Ambient	

Formulation

Component	concentration mg/l
D-Glucose	1000
KCI	400
KH ₂ PO ₄	60
NaCl	8000
NaHCO ₃	350
Na ₂ HPO ₄	48
Phenol red sodium salt	11

Precautions and Disclaimer

Thisproductisforresearchuseonly.