

## Datasheet

<b>Catalogue #</b>	LT-AG-SA-01
<b>Product name</b>	Streptavidin recombinant
<b>Source</b>	from recombinant E. coli strain
<b>Specificity</b>	React specifically with biotin containing biomolecules
<b>Molecular weight</b>	19,2 kD
<b>Theoretical pI</b>	6,43
<b>Purification tag</b>	6xHis
<b>Unit</b>	1mg, 10 mg
<b>Concentration</b>	Lot dependent
<b>Purification</b>	chromatography
<b>Buffer</b>	PBS
<b>Stabilizer</b>	none, may be added on request
<b>Preservative</b>	None
<b>Application</b>	ELISA, Lateral Flow
<b>Sterility</b>	Filter sterilized
<b>Long Term Storage</b>	Lyophilized
<b>Material Safety note</b>	This product is sold as an antigen preparation for research use only. Standard Laboratory Practices should be followed when handling this material.
<b>Validation data:</b>	

SDS-PAGE

